




Freedom™
Trans-Obturator Sling



Freedom™
Trans-Obturator Needles

Freedom™


Trans-Obturator Needles and Sling System

for Female Stress Urinary Incontinence


Product Description:

The Freedom™ Trans Obturator Sling System is intended for the placement of a sub-urethral mesh for the treatment of female stress urinary incontinence (SUI) resulting from urethral hyper mobility and/or intrinsic sphincter deficiency.

Obturator Anatomy:



Product Details:

Sling:	Monofilament Dia. mm (in)	Typical Pore Size mm (in)	Burst Strength kpa (psi)	Break Strength N/2.5 cm (lbs/in)		Break Elongation (%)		Thick mm (in)	Wt. g/m ² (oz/yd ²)
				MD	CMD	MD	CMD		
	0.15 (.006)	1.3 x 1.0 (.051 x .039)	1172 (162)	526 (118)	492 (111)	107	175	0.59 (.023)	100 (2.9)

Code	Description
F VM-S	Freedom Trans-Obturator Sling System
F RM-NDL	1 Set of Freedom Trans-Obturator Needles with CD
F RM-KIT	Kit of 1 Set Freedom Needles + 3 Freedom Trans-Obturator Sling System + CD